

Table of Contents

data download and setup.....	1
load metabolite structure	1
Inputs and flags for metabolite mapping.....	1
Process data.....	1
Plotting.....	1
Data structure of Metcon_CSI object.....	1
Processing CSI-bSSFP data.....	2
Other miscellaneous methods which can be useful.....	2
Debug (only for AMARES and Lorentzian solver so far).....	2

Demo 01: 3D CSI Phantom data

data download and setup

add paths and set up data paths. Download phantom data (1.82 GB "phantom-DMI.tar.gz") from [zenodo.org](#) (DOI: 10.5281/zenodo.14652737).

```
% add all dependencies
addpath(genpath('/ptmp/pvalsala/Packages/mapVBVD'))
addpath(genpath('/ptmp/pvalsala/Packages/DeuteMetCon'))
addpath(genpath('/ptmp/pvalsala/Packages/OXSA'))

% data path
sn='/ptmp/pvalsala/deuterium/dataForPublication/phantom-DMI';
dirst_csi=dir(fullfile(sn,"*rpcsi_fid*.dat"));
```

load metabolite structure

The chemical shifts, measured relaxation times and labels of all four metabolites in phantom is organised in a array of struct

```
% second argument is the frequency offset of water in Hz
metabolites=getMetaboliteStruct('phantom',0);
```

Inputs and flags for metabolite mapping

All functions and data required for all data processing steps like image reconstruciton, coil combination and spectral seperation were encapsualted in a single class `MetCon_CSI.m`. All inputs and flags except the raw data file is a name-value pair as described below.

%	name	description	default	possible options
%	-----	-----	-----	-----
%	'metabolites'	struct array with definition of metabolites	[]	see getMetaboliteStruct.m function
%	'fm'	1H fieldmap in rad/s	[]	3D numeric matrix or 'IDEAL'
%	'csm'	coil maps	[]	3D numeric matrix
%	'mask'	mask for spectral separation	[]	3D logical matrix
%	'doDenosing'	SVD denoising	0	scalar No of components, -1 for debug
%	'Solver'	spectral separation method	'IDEAL'	{'phaseonly','pinv','IDEAL','IDEAL-modes','AMARES','LorentzFit'}
%				'phaseonly'- linear method with only phase evolution
%				'pinv'- linear method with full signal model
%				'IDEAL'- iterative IDEAL algorithm
%				'IDEAL-modes'-IDEAL algorithm for phase cycled data
%				'AMARES'- AMARES spectral fitting
%				'LorentzFit'- lorentzian spectral fitting
%	'parfor'	flag to use parfor	true	boolean
%	'doZeroPad'	zero pad factor	[1 1 1 0]	positive scalar array [3 physical axis x 1 time]
%	'doSmoothFM','maxit'	IDEAL flags: fieldmap smooth factor and maximum iterations	1,10	scalar(+ve: gaussian, -ve: median),postive scalar
%	'doPhaseCorr'	phase correction mode	'none'	{'none','Manual','Burg'}
%	'CoilSel','PCSel','EchoSel'	arrays to picks some of coils, time points and phasecycles.	1:max()	positive integer array
%	'doNoiseDecorr'	flag to perform noise decorrelation	true	boolean
%	'doCoilCombine'	coil combine mode	'adapt1'	{'none','sos','adapt1','wsvd'}

```
CSI_setting={'metabolites',metabolites,'doPhaseCorr','none','parfor',true,...
            'doCoilCombine','adapt1','doZeroPad',[0.5 0.5 0.5 0],'mask',[],'Solver','IDEAL','fm',[]};
```

Process data


```
CSI_filename=fullfile(sn,dirst_csi(end).name);
mcobj_csi=MetCon_CSI(CSI_filename,CSI_setting{:});
```

Software version: VD (1?)
Reader version: 1660732089 (UTC: 17-Aug-2022 10:28:09)
Scan 1/1, read all mdhs:
 661.3 MB read in 10 s
Using T2* =20 ms
starting reco
initial CSI data size: 10 25 25 25 256
final CSI data size: 10 51 51 51 256
reco time = 19.8 s
estimating field map(1/2)
estimating metabolities(2/2)
Metabolite mapping time = 26.8 s

Plotting

After image reconstruciton and spectral seperation, metbolite amplitudes are store in mcobj_csi.Metcon. The 4D Metcon matrix and 2H field map estimated by IDEAL algorithm can be quickly visualized with `PlotResults` method.

```
mcobj_csi.PlotResults()
```



Data structure of Metcon_CSI object

Property	Value
time	[]
FieldMap	[]
mask	51x51x51 logical
twix	1x1 struct
DMIPara	1x1 struct
flags	1x1 struct
filename	/ptmp/pvalsala/deuterium/dataForPublication/phantom/MID01423...
metabolites	1x4 struct
sig	5-D complex single
img	5-D complex single
coilSens	4-D complex single
coilNormMat	51x51x51 single
SolverObj	1x1 IDEAL
Metcon	4-D complex double
Experimental	1x1 struct
D	10x10 complex single

Where

- DMIPara are the important sequence parameters parsed from twix with getDMIPara.m fucntion
- flags contains all the processing flags
- sig - averaged signal [CHAx LIN xSEG x PAR x COL (time) x REP(PhaseCylces)]
- img - reconstructed image [CHA x 3 physical dimension x time]
- Metcon - Metabolite amplitudes [CHA x 3 physical dimension x time]
- Experimental- contains all experimental outputs of processing (fieldmap, residue, other fit parameters, fit quality)
- D- noise decorrelation matrix.

Processing CSI-bSSFP data


We use `IDEAL-modes` for phase-cycled data for SNR optimallity. Other Solver flags use averaged phase-cyles volumes

```
CSI_setting_ssfp={ 'metabolites',metabolites, 'doPhaseCorr', 'none', 'parfor',true,...
    'doCoilCombine', 'adapt1', 'doZeropad',[0.5 0.5 0.5 0], 'mask',[], 'Solver', 'IDEAL-modes', 'fm', 'IDEAL'};
dirst_csi_ssfp=dir(fullfile(sn, "**pcsi_ssfp*.dat"));
CSI_filename=fullfile(sn,dirst_csi_ssfp(end).name);
mcobj_ssfp=MetCon_CSI(CSI_filename,CSI_setting_ssfp{:});

Software version: VD (!?)
Reader version: 1660732089 (UTC: 17-Aug-2022 10:28:09)
Scan 1/1, read all mdhs:
    335.9 MB read in   10 s
Using T2* =20 ms
starting reco
initial CSI data size:      10 25 25 25 64 4
final CSI data size:      10 51 51 51 64 4
reco time =   17.7 s
estimating field map(1/2)
estimating metabolities(2/2)
Metabolite mapping time =   25.1 s

mcobj_ssfp.PlotResults()
```

M1425|TR 19 ms| 48 deg | 8.32 mm | 4 rep | 64 echoes|IDEAL-mo




Other miscallaneous methods which can be useful

```
% The metabolite amplitudes can be normalized into SNR
metcon_SNR= mcobj_csi.getNormalized();

% method and quantified into mM with 10 mM water reference using `getmM`
% fuction. (Doesn't work in phantom as 10mM water reference is not there!)
metcon_mM= mcobj_csi.getmM();

% export NIFTI volumes of all outputs (average image, metabolite amplitude in SNR unit and mM)
pn=fullfile(sn, 'proc',sprintf('CSI_%s',datetime('today','Format','yyyyMMdd')));
mkdir(pn);cd(pn);
niiFileName=mcobj_csi.WriteImages(pn);

%plot the k-space weighting and PSF to get realistic voxel size (FW64%)
voxel_size_mm= getPSF_CSI(mcobj_csi.twix);
```



```
% it is non-sense for phantom scans but still!
Intake_time_mins=mcobj_csi.getMinutesAfterIntake('08:00');
```


Debug (only for AMARES and Lorentzian solver so far)

Useful for checking data at a particular voxel index. [15,14,25] is the voxel with all metabolites.

```
mcobj_csi.demoFit([15,14,25])

Iterations = 36.
Norm of residual = 770.811
Norm of the data = 62474.292
resNormSq / dataNormSq = 0.012
zeroOrderPhaseRad = 2.222
    chemShift: [0.0474 -1.0161 -2.4680 -3.4884]
    linewidth: [2.0000 11.2813 4.7048 2.2230]
    amplitude: [34.5440 11.3873 8.7223 14.5698]
           phase: [127.3026 127.3026 127.3026 127.3026]
           sigma: [0 0 0 0]

Performing Nlorentz fit
```



Warning: Converting Y to vector of double.
Fitting stopped because the number of iterations or function evaluations exceeded the specified maximum.
Warning: Imaginary parts of complex X and/or Y arguments ignored.

