

Table of Contents	
data download and setup.....	1
load metabolite structure	1
Inputs and flags for metabolite mapping.....	1
Process data.....	1
Plotting.....	1
Data structure of Metcon_CSI object.....	1
Processing CSI-bSSFP data.....	2
Overlay plot with anatomy.....	2
Other miscallaneous methods which can be useful.....	2
Debug (only for AMARES and Lorentzian solver so far).....	3

Demo 01: 3D CSI Phantom data

data download and setup

add paths and set up data paths. Download in vivo data (3.8 GB "sub-01-DMI.tar.gz") from [zenodo.org](https://zenodo.org/record/14652737/files/sub-01-DMI.tar.gz) (DOI: 10.5281/zenodo.14652737).

```
% add all dependencies
addpath(genpath('/ptmp/pvalsala/Packages/mapVBVD'))
addpath(genpath('/ptmp/pvalsala/Packages/DeuteMetCon'))
addpath(genpath('/ptmp/pvalsala/Packages/OXSA'))

% set up data folder
sn='/ptmp/pvalsala/deuterium/dataForPublication/sub-01-DMI';
dirst_csi=dir(fullfile(sn,"*rpcsi_fid*.dat"));
```

load metabolite structure

The chemical shifts, measured relaxation times and labels of all four metabolites in phantom is organised in a array of struct

```
% second argument is the frequency offset of water in Hz
metabolites=getMetaboliteStruct('invivo',0);
```

Inputs and flags for metabolite mapping

All functions and data required for all data processing steps like image reconstruciton, coil combination and spectral seperation were encapsualted in a single class `MetCon_CSI.m`. All inputs and flags except the raw data file is a name-value pair as described below.

%	name	description	default	possible options
%	-----	-----	-----	-----
%	'metabolites'	struct array with definition of metabolites	[]	see getMetaboliteStruct.m function
%	'fm'	1H fieldmap in rad/s	[]	3D numeric matrix or 'IDEAL'
%	'csm'	coil maps	[]	3D numeric matrix
%	'mask'	mask for spectral separation	[]	3D logical matrix
%	'doDenosing'	SVD denoising	0	scalar No of components, -1 for debug
%	'Solver'	spectral separation method	'IDEAL'	{'phaseonly','pinv','IDEAL','IDEAL-modes','AMARES','LorentzFit'}
%				'phaseonly'- linear method with only phase evolution
%				'pinv'- linear method with full signal model
%				'IDEAL'- iterative IDEAL algorithm
%				'IDEAL-modes'-IDEAL algorithm for phase cycled data
%				'AMARES'- AMARES spectral fitting
%				'LorentzFit'- lorentzian spectral fitting
%				
%	'parfor'	flag to use parfor	true	boolean
%	'doZeroPad'	zero pad factor	[1 1 1 0]	positive scalar array [3 physical axis x 1 time]
%	'doSmoothFM','maxit'	IDEAL flags: fieldmap smooth factor and maximum iterations	1,10	scalar(+ve: gaussian, -ve: median),postive scalar
%	'doPhaseCorr'	phase correction mode	'none'	{'none','Manual','Burg'}
%	'CoilSel','PCSel','EchoSel'	arrays to picks some of coils, time points and phasecycles.	1:max()	positive integer array
%	'doNoiseDecorr'	flag to perform noise decorrelation	true	boolean
%	'doCoilCombine'	coil combine mode	'adapt1'	{'none','sos','adapt1','wsvd'}

```
CSI_setting={'metabolites',metabolites,'doPhaseCorr','none','parfor',true,...
            'doCoilCombine','adapt1','doZeropad',[0.5 0.5 0.5 0],'mask',[],'Solver','IDEAL','fm',[]};
```

Process data

```
CSI_filename=fullfile(sn,dirst_csi(end).name);
mcobj_fisp=MetCon_CSI(CSI_filename,CSI_setting{:});
```

```
Software version: VD (!?)
Reader version: 1660732089 (UTC: 17-Aug-2022 10:28:09)
Scan 1/1, read all mdhs:
    661.3 MB read in    8 s
Using T2* =20 ms
starting reco
initial CSI data size:      10   25   25   25  256
final CSI data size:      10   51   51   51  256
reco time =   20.1 s
estimating field map(1/2)
estimating metabolities(2/2)
Metabolite mapping time =   28.2 s
```

Ploting

After image reconstruciton and spectral seperation, metbolite amplitudes are store in mcobj_csi.Metcon. The 4D Metcon matrix and 2H field map estimated by IDEAL algorithm can be quickly visualized with `PlotResults` method.

```
mcobj_fisp.PlotResults()
```



Data structure of Metcon_CSI object

1x1 MetCon_CSI	
Property	Value
time	[]
FieldMap	[]
mask	51x51x51 logical
twix	1x1 struct
DMIPara	1x1 struct
flags	1x1 struct
filename	'/ptmp/pvalsala/deuterium/dataForPublication/phantom/MID01423...
metabolites	1x4 struct
sig	5-D complex single
img	5-D complex single
coilSens	4-D complex single
coilNormMat	51x51x51 single
SolverObj	1x1 IDEAL
Metcon	4-D complex double
Experimental	1x1 struct
D	10x10 complex single

Where

- DMIPara are the important sequence parameters parsed from twix with getDMIPara.m fucntion
- flags contains all the processing flags
- sig - averaged signal [CHAx LIN xSEG x PAR x COL (time) x REP(PhaseCylces)]
- img - reconstructed image [CHA x 3 physical dimension x time]
- Metcon - Metabolite amplitudes [CHA x 3 physical dimension x time]
- Experimental- contains all experimental outputs of processing (fieldmap, residue, other fit parameters, fit quality)
- D- noise decorrelation matrix.

Processing CSI-bSSFP data


We use `IDEAL-modes` for phase-cycled data for SNR optimallity. Other Solver flags use averaged phase-cyles volumes

```
CSI_setting_ssfp={'metabolites',metabolites,'doPhaseCorr','none','parfor',true,...
    'doCoilCombine','adapt1','doZeropad',[0.5 0.5 0.5 0],'mask',[],'Solver','IDEAL-modes','fm','IDEAL'};
dirst_csi_ssfp=dir(fullfile(sn,"*rpcsi_ssfp*.dat"));
CSI_filename=fullfile(sn,dirst_csi_ssfp(1).name);
mcobj_ssfp=MetCon_CSI(CSI_filename,CSI_setting_ssfp{:});

Software version: VD (!?)
Reader version: 1660732089 (UTC: 17-Aug-2022 10:28:09)
Scan 1/1, read all mdhs:
    335.9 MB read in    7 s
Using T2* =20 ms
starting reco
initial CSI data size:      10 25 25 25 64 4
final CSI data size:      10 51 51 51 64 4
reco time = 17.4 s
estimating field map(1/2)
estimating metabolities(2/2)
Metabolite mapping time = 25.5 s

mcobj_ssfp.PlotResults()
```

M1000|TR 19 ms| 48 deg | 8.32 mm | 4 rep | 64 echoes|IDEAL-mo



Overlay plot with anatomy

we need spm12 in path for reslicing. We plot the third metabolite (Glx) over skull stripped anatomy.

```
%addpath('/ptmp/pvalsala/Packages/spm12/')
anat_tra=fullfile(sn,'anat/brain_tra.nii');

% export NIFTI volumes of all outputs (average image, metabolite amplitude in SNR unit and mM)
pn=fullfile(sn,sprintf('proc/csi_%s',datetime('today','Format','yyyyMMdd')));
mkdir(pn);
niiFileName_fisp=mcobj_fisp.WriteImages(pn);
niiFileName_ssfp=mcobj_ssfp.WriteImages(pn);
% reslice the metabolite amplitude to anatomy!
resliced_metcon=myspm_reslice(anat_tra,dir(fullfile(pn,'Metcon_SNR_*csi*.nii')),'nearest','rt');
```


```
Warning: Trying to add spm12
reslicing 2 volumes: Metcon_SNR_m00997_rpcsi_fid_Stan_res15.6_moreoptimal_IDEAL_none.nii
Metcon_SNR_m01000_rpcsi_ssfp_Stan25_15.6mm_IDEAL-modes_none.nii

-----
23-Jan-2025 18:32:05 - Running job #1
-----
23-Jan-2025 18:32:05 - Running 'Coregister: Reslice'

SPM12: spm_reslice (v7141) 18:32:06 - 23/01/2025
=====
Completed : 18:32:08 - 23/01/2025
23-Jan-2025 18:32:08 - Done 'Coregister: Reslice'
23-Jan-2025 18:32:08 - Done

% Overlay plot
clim_glx=[0 30];
figure,
overlayplot(dir(anat_tra),dir(fullfile(pn,'rt*.nii')),'MetIdx',3,'SlcSel',15, ...
    'transform',@(x) flippermute(x,[2 1 3 4 5]),1),...
    'cax',clim_glx,'cmap','turbo');
title([metabolites(3).name,' [SNR]'] % 'MetIdx',3
y ticklabels({sprintf('FISP %d mins',mcobj_fisp.getMinutesAfterIntake('08:36')),...
    sprintf('bSSFP %d mins',mcobj_ssfp.getMinutesAfterIntake('08:36'))});

%make your own colorbar
cb_handle=colorbar;
cb_handle.Visible='off';
ax2 = axes('Position',cb_handle.Position);
imagesc(linspace(0,1,100)'),colormap(ax2,'turbo')
set(ax2,'YAxisLocation','right','FontSize',10,'FontWeight','bold','YDir','normal')
yticks(round(linspace(0,100,4)))
yticklabels(round(linspace(0,1,4)*clim_glx(2),1))
xticks([])
```



Other miscallaneous methods which can be useful


```
% The metabolite amplitudes can be normalized into SNR
metcon_SNR= mcobj_fisp.getNormalized();

% method and quantified into mM with 10 mM water reference using `getmM`
% fucntion. (Doesn't work in phantom as 10mM water reference is not there!)
metcon_mM= mcobj_fisp.getmM();
```

```
% export NIFTI volumes of all outputs (average image, metabolite amplitude in SNR unit and mM)
pn=fullfile(sn,sprintf('proc/csi_%s',datetime('today','Format','yyyyMMdd')));
niiFileName=mcobj_fisp.WriteImages(pn);

%plot the k-space weighting and PSF to get realistic voxel size (FW64%)
voxel_size_mm= getPSF_CSI(mcobj_fisp.twix);
```

Using T2* =20 ms




```
% it is non-sense for phantom scans but still!
Intake_time_mins=mcobj_fisp.getMinutesAfterIntake('08:00');
```

Debug (only for AMARES and Lorentzian solver so far)

Useful for checking data at a particular voxel index. [38,30,34] is gray matter voxel.

```
mcobj_fisp.demoFit([38,30,34])
```



Performing Nlorentz fit
Warning: Converting Y to vector of double.
Fitting stopped because the number of iterations or function evaluations exceeded the specified maximum.
Warning: Imaginary parts of complex X and/or Y arguments ignored.

