

Table of Contents

data download and setup.....	1
load metabolite structure	1
Inputs and flags for metabolite mapping.....	1
Get noise decoorelation matrix	1
assemble inputs and processing flags.....	1
Process data.....	1
Plotting.....	1
Data structure of Metcon_ME object.....	1
Other miscallaneous methods which can be useful.....	2
Debug (only for 'pinv' (linear) fit mode).....	2

Demo 02: 3D ME Phantom data

data download and setup

add paths and set up data paths. Download phantom data (1.82 GB "phantom-DML.tar.gz") from [zenodo.org](#) (DOI: 10.5281/zenodo.14652737).

```
% add all dependencies
addpath(genpath('/ptmp/pvalsala/Packages/mapVBVD'))
addpath(genpath('/ptmp/pvalsala/Packages/DeuteMetCon'))

% data path
sn='/ptmp/pvalsala/deuterium/dataForPublication/phantom-DMI';
dirst_me=dir(fullfile(sn,"*trufi_5E*.dat"));
dirst_noise=dir(fullfile(sn,"*trufi*noise*.dat"));
```

load metabolite structure

The chemical shifts, measured relaxation times and labels of all four metabolites in phantom is organised in a array of struct

```
% second argument is the frequency offset of water in Hz
metabolites=getMetaboliteStruct('phantom',0);
```

Inputs and flags for metabolite mapping

All functions and data required for all data processing steps like image reconstruciton, coil combination and spectral seperation were encapsualted in a single class `MetCon_CSI.m`. All inputs and flags except the raw data file is a name-value pair as described below.

% name	description	default	possible options	
% -----	-----	-----	-----	
% metabolites	struct array with definition of metabolites	[]	see getMetaboliteStruct.m function	
% fm	1H fieldmap in rad/s	[]	3D numeric matrix or 'IDEAL'	
% csm	coil maps	[]	3D numeric matrix	
% mask	mask for spectral separation	[]	3D logical matrix , scalar percentile (1-100) threshold	
% doDenosing	SVD denoising	0	scalar No of components, -1 for debug	
% Solver	spectral separation method	'IDEAL'	{'phaseonly','pinv','IDEAL','IDEAL-modes','AMARES','LorentzFit'}	
%			'phaseonly'- linear method with only phase evolution	
%			'pinv'- linear method with full signal model	
%			'IDEAL'- iterative IDEAL algorithm	
%			'IDEAL-modes'-IDEAL algorithm for phase cycled data	
%			'AMARES'- AMARES spectral fitting	
%			'LorentzFit'- lorentzian spectral fitting	
% parfor	flag to use parfor	true	boolean	
% doZeroPad	zero pad factor	[1 1 1 0]	positive scalar array [3 physical axis x 1 time]	
% 'doSmoothFM','maxit'	IDEAL flags: fieldmap smooth factor and maximum iterations	1,10	scalar(+ve: gaussian, -ve: median),postive scalar	
% doPhaseCorr	phase correction mode	'none'	{'none','Manual','Burg'}	
% 'CoilSel','PCSel','EchoSel'	arrays to picks some of coils, time points and phasecycles.	1:max()	positive integer array	
% doNoiseDecorr	flag to perform noise decorrelation	true	boolean	

Get noise decoorelation matrix

trufi sequence acquire noise data only when parallel imaging is enabled. Therefore, we acquire noise scan with 0 flip angle seperately.

```
twix_noise=mapVBVD(fullfile(sn,dirst_noise(1).name),'rmos');

Software version: VD (1?)
Reader version: 1660732089 (UTC: 17-Aug-2022 10:28:09)
Scan 1/1, read all mdhs:
    40.4 MB read in    6 s

[D_noise,D_image,noise_info]=CalcNoiseDecorrMat(twix_noise);
%D_image is the noise correlation from image data measured with 0 FA
```

assemble inputs and processing flags

```
ME_setting={ 'NoiseDecorr',D_image,'mask',[],'metabolites',metabolites,...
             'doPhaseCorr',true,'doZeropad',[1 1 1]*0.5,'parfor',true,'fm','IDEAL','Solver','pinv'};

%
%'fm','IDEAL' : estimate field map from the averaged phasecycle volume using IDEAL algorithm
%'Solver','pinv' : linear fit (least square fit of full signal model)
```

Process data

```
ME_filename=fullfile(sn,dirst_me(end).name);
mcobj_me=MetCon_ME(ME_filename,ME_setting{:});

Software version: VD (1?)
Reader version: 1660732089 (UTC: 17-Aug-2022 10:28:09)
Scan 1/1, read all mdhs:
    856.4 MB read in    8 s
starting reco
estimating field map(1/2)
estimating metabolities(2/2)
reco time =   19.0 s
Calculating bSSFP profile basis
done.....


Metabolite fitting done in 5.2 s !
Performed pixel shift along read: (0.0,-2.2,-5.2 , -7.2) mm
Metabolite mapping time =    5.3 s
```

Plotting

















After image reconstruciton and spectral seperation, metbolite amplitudes are store in mcobj_csi.Metcon. The 4D Metcon matrix and 2H field map estimated by IDEAL algorithm can be quickly visualized with `PlotResults` method.

```
mcobj_me.PlotResults()
```

M1428|TR 19 ms| 48 deg | 12.50 mm | 18 rep | 5 echoes|pinv



Data structure of Metcon_ME object

 1x1 MetCon_ME	
Property	Value
 FieldMap	48x64x32 double
 mask	48x64x32 logical
 twix	1x1 struct
 DMIPara	1x1 struct
 flags	1x1 struct
 filename	/ptmp/pvalsala/deuterium/dataForPublication/phantom-DM/MID0...
 metabolites	1x4 struct
 sig	6-D complex single
 img	6-D complex single
 coilSens	4-D complex single
 coilNormMat	48x64x32 single
 SolverObj	[]
 Metcon	4-D complex double
 Experimental	1x1 struct
 D	10x10 complex single

Where

- DMIPara are the important sequence parameters parsed from twix with `getDMIPara.m` fucntion
- flags contains all the processing flags
- sig - averaged signal [CHA x LIN x COL x PAR x ECO x REP]
- img - reconstructed image [CHA x Phase x Read xSlice x echo x PC]
- Metcon - Metabolite amplitudes [CHA x 3 physical dimension x time]
- Experimental- contains all experimental outputs of processing (fieldmap, residue, other fit parameters, fit quality)
- D- noise decorrelation matrix.

Other miscallaneous methods which can be useful


```
% The metabolite amplitudes can be normalized into SNR
metcon_SNR= mcobj_me.getNormalized();

% method and quantified into mM with 10 mM water reference using `getmM`
% fuction.
metcon_mM= mcobj_me.getMm();

% caculate measurement time after glucose intake in mins
Intake_time_mins=mcobj_me.getMinutesAfterIntake('08:00');

% export NIFTI volumes of all outputs (average image, metabolite amplitude in SNR unit and mM)
pn=fullfile(sn,'proc',sprintf('ME_%s',datetime('today','Format','yyyyMMdd')));
mkdir(pn);
niiFileName=mcobj_me.WriteImages(pn);

%plot the k-space weighting and PSF to get realistic voxel size (FW64%)
voxel_size_mm= getPSF_CSI(mcobj_me.twix,true)
```



```
voxel_size_mm = 4x1
    12.4793
    14.4160
    14.8232
         0
```

Debug (only for 'pinv' (linear) fit mode)

Useful for checking data at a particular voxel index. [24,33,27] is a lactate voxel.

```
%prints metabolite amplitudes and plot fit and basis functions
mcobj_me.demoFit([24,33,27])
```

```
12.4828
1.6200
0.1262
77.7509
```

