

Table of Contents

data download and setup.....	1
load metabolite structure .....	1
Inputs and flags for metabolite mapping.....	1
Get noise decoorelation matrix .....	1
assemble inputs and processing flags.....	1
Process data.....	1
Plotting.....	1
Data structure of Metcon_ME object.....	1
Other miscallaneous methods which can be useful.....	2
Debug (only for 'pinv' (linear) fit mode).....	2

Demo 02: 3D ME Phantom data

data download and setup

add paths and set up data paths. Download phantom data (1.82 GB "phantom-DML.tar.gz") from [zenodo.org](https://zenodo.org/record/14652737/files/phantom-DML.tar.gz) (DOI: 10.5281/zenodo.14652737).

```
% add all dependencies
addpath(genpath('/ptmp/pvalsala/Packages/mapVBVD'))
addpath(genpath('/ptmp/pvalsala/Packages/DeuteMetCon'))

% data path
sn='/ptmp/pvalsala/deuterium/dataForPublication/phantom-DMI';
dirst_me=dir(fullfile(sn,"*trufi_5E*.dat"));
dirst_noise=dir(fullfile(sn,"*trufi*noise*.dat"));
```

load metabolite structure

The chemical shifts, measured relaxation times and labels of all four metabolites in phantom is organised in a array of struct

```
% second argument is the frequency offset of water in Hz
metabolites=getMetaboliteStruct('phantom',0);
```

Inputs and flags for metabolite mapping

All functions and data required for all data processing steps like image reconstruciton, coil combination and spectral seperation were encapsualted in a single class `MetCon\_CSI.m`. All inputs and flags except the raw data file is a name-value pair as described below.

%   name	description	default	possible options	
%   -----	-----	-----	-----	
%   metabolites	struct array with definition of metabolites	[ ]	see getMetaboliteStruct.m function	
%   fm	1H fieldmap in rad/s	[ ]	3D numeric matrix or 'IDEAL'	
%   csm	coil maps	[ ]	3D numeric matrix	
%   mask	mask for spectral separation	[ ]	3D logical matrix , scalar percentile (1-100) threshold	
%   doDenosing	SVD denoising	0	scalar No of components, -1 for debug	
%   Solver	spectral separation method	'IDEAL'	{'phaseonly','pinv','IDEAL','IDEAL-modes','AMARES','LorentzFit'}	
%			'phaseonly'- linear method with only phase evolution	
%			'pinv'- linear method with full signal model	
%			'IDEAL'- iterative IDEAL algorithm	
%			'IDEAL-modes'-IDEAL algorithm for phase cycled data	
%			'AMARES'- AMARES spectral fitting	
%			'LorentzFit'- lorentzian spectral fitting	
%   parfor	flag to use parfor	true	boolean	
%   doZeroPad	zero pad factor	[1 1 1 0]	positive scalar array [3 physical axis x 1 time]	
%   'doSmoothFM','maxit'	IDEAL flags: fieldmap smooth factor and maximum iterations	1,10	scalar(+ve: gaussian, -ve: median),postive scalar	
%   doPhaseCorr	phase correction mode	'none'	{'none','Manual','Burg'}	
%   'CoilSel','PCSel','EchoSel'	arrays to picks some of coils, time points and phasecycles.	1:max()	positive integer array	
%   doNoiseDecorr	flag to perform noise decorrelation	true	boolean	

Get noise decoorelation matrix

trufi sequence acquire noise data only when parallel imaging is enabled. Therefore, we acquire noise scan with 0 flip angle seperately.

```
twix_noise=mapVBVD(fullfile(sn,dirst_noise(1).name),'rmos');
```

```
Software version: VD (1?)
Reader version: 1660732089 (UTC: 17-Aug-2022 10:28:09)
Scan 1/1, read all mdhs:
Warning: MATLAB has disabled some advanced graphics rendering features by switching to software OpenGL. For more information, click here.
40.4 MB read in 8 s
```

```
[D_noise,D_image,noise_info]=CalcNoiseDecorrMat(twix_noise);
%D_image is the noise correlation from image data measured with 0 FA
```

assemble inputs and processing flags

```
ME_setting={ 'NoiseDecorr',D_image,'mask',[ ],'metabolites',metabolites,...
             'doPhaseCorr',true,'doZeropad',[1 1 1]*0.5,'parfor',true,'fm','IDEAL','Solver','pinv'};
```

```
%
%'fm','IDEAL' : estimate field map from the averaged phasecycle volume using IDEAL algorithm
%'Solver','pinv' : linear fit (least square fit of full signal model)
```

Process data


```
ME_filename=fullfile(sn,dirst_me(end).name);
mcobj_me=MetCon_ME(ME_filename,ME_setting{:});
```

```
Software version: VD (1?)
Reader version: 1660732089 (UTC: 17-Aug-2022 10:28:09)
Scan 1/1, read all mdhs:
856.4 MB read in 8 s
starting reco
estimating field map(1/2)
Starting parallel pool (parpool) using the 'Processes' profile ...
Warning: The system time zone setting, 'CET', does not specify a single time zone unambiguously. It will be treated as 'Europe/Zurich'. See the datetime.TimeZone property for details about specifying time zones.
Connected to parallel pool with 36 workers.
estimating metabolities(2/2)
reco time = 65.2 s
Calculating bSSFP profile basis
done.....

Metabolite fitting done in 5.2 s !
Performed pixel shift along read: (0.0,-2.2,-5.2 , -7.2) mm
Metabolite mapping time = 5.3 s
```

Plotting

After image reconstruciton and spectral seperation, metbolite amplitudes are store in mcobj\_csi.Metcon. The 4D Metcon matrix and 2H field map estimated by IDEAL algorithm can be quickly visualized with `PlotResults` method.



Data structure of Metcon\_ME object

1x1 [MetCon\\_ME](#)

Property	Value
FieldMap	48x64x32 double
mask	48x64x32 logical
twix	1x1 struct
DMIPara	1x1 struct
flags	1x1 struct
filename	/ptmp/pvalsala/deuterium/dataForPublication/phantom-DM/MID0...
metabolites	1x4 struct
sig	6-D complex single
img	6-D complex single
coilSens	4-D complex single
coilNormMat	48x64x32 single
SolverObj	[ ]
Metcon	4-D complex double
Experimental	1x1 struct
D	10x10 complex single

Where

- DMIPara are the important sequence parameters parsed from twix with `getDMIPara.m` fucntion
- flags contains all the processing flags
- sig - averaged signal [CHA x LIN x COL x PAR x ECO x REP]
- img - reconstructed image [CHA x Phase x Read xSlice x echo x PC]
- Metcon - Metabolite amplitudes [CHA x 3 physical dimension x time]
- Experimental- contains all experimental outputs of processing (fieldmap, residue, other fit parameters, fit quality)
- D- noise decorrelation matrix.

Other miscallaneous methods which can be useful


```
% The metabolite amplitudes can be normalized into SNR
metcon_SNR= mcobj_me.getNormalized();

% method and quantified into mM with 10 mM water reference using `getmM`
% fuction.
metcon_mM= mcobj_me.getM();

% caculate measurement time after glucose intake in mins
Intake_time_mins=mcobj_me.getMinutesAfterIntake('08:00');

% export NIFTI volumes of all outputs (average image, metabolite amplitude in SNR unit and mM)
pn=fullfile(sn,'proc',sprintf('ME_%s',datetime('today','Format','yyyyMMdd')));
mkdir(pn);
niiFileName=mcobj_me.WriteImages(pn);

%plot the k-space weighting and PSF to get realistic voxel size (FW64%)
voxel_size_mm= getPSF_CSI(mcobj_me.twix)
```



Debug (only for 'pinv' (linear) fit mode)

Useful for checking data at a particular voxel index. [24,33,27] is a lactate voxel.

```
%prints metabolite amplitudes and plot fit and basis functions
mcobj_me.demoFit([24,33,27])
```

12.4828  
1.6200  
0.1262  
77.7509

